

	<b>Technical data sheet</b>	Ref : SBM1 Page : 1/2 Version date : 04/06/19
---	-----------------------------	---

## OPTIMUS SMG TAQ DNA POLYMERASE

A heat-stable DNA polymerase with a molecular weight of 95kDa suitable for routine PCR applications, standard and high throughput testing, and screening. It has 5' to 3' exonuclease activity to leave 3' end overhang with adenine bases for direct cloning. Available in 500 units per transaction.

**CAT No.:** SBM1

**Molecular Weight:** 95kDa

**Concentration:** 5 units/ $\mu$ l

**Application:**

- Routine PCR
- Colony PCR
- Standard testing
- Screening
- High throughput testing

**Transportation:** Kept on blue ice

**Content:** Storage Buffer

- 25mM Tris-HCl (pH8.0), 100mM KCl, 0.1mM EDTA, 1mM DTT, 50% Glycerol, 0.5% Nonident P40, 0.5% Tween 20

**Reaction Buffers supplied with the enzyme:**

- 10X Buffer I: 500mM KCl, 100mM Tris-HCl, pH 9.0, 1% Triton X-100, 15mM MgCl<sub>2</sub>
- 10X Buffer II: 500mM KCl, 100mM Tris-HCl, p H 9.0. 1% Triton X-100
- MgCl<sub>2</sub>: 25 mM

**Usage:**

- Components Volume per reaction:

- 10X reaction buffer I or buffer II 5  $\mu$ l
- 25 mM MgCl<sub>2</sub> 1.5  $\mu$ l (if you use buffer II)
- dNTP-Mix (40mM) 1.0  $\mu$ l
- Up-stream primer (10  $\mu$ M stock) 0,5-2.5  $\mu$ l

- Down-stream primer (10 $\mu$ M stock) 0.5-2,5  $\mu$ l
- Template DNA  
0.1-15 ng/ml plasmid DNA
- 1-10  $\mu$ g/ml genomic DNA
- Maximo Taq DNA Polymerase (5 u/ $\mu$ l) 0.2 - 1.0  $\mu$ l
- Sterile dest. Water (molecular grade) up to 50  $\mu$ l total reaction volume

**Storage:** Store at -20°C for 24 months

**Quality Control:**

- PCR with various templates—genomic DNA, Phage Lambda DNA
- 3kb DNA amplification from 50ng DNA
- Batch variation and level of bacterial DNA contamination