

	Technical data sheet	Ref : SMGBM-1001RL
		Page : 1/1
		Version date : 04/06/19

SMG MAS DNA MARKER 100bp

SMG MAS DNA marker provides the most accurate sizing of fragments over a broad range of sizes. Can be loaded into agarose gels. Simplifies and estimates the molecular weight of double-stranded DNA fragments generated by PCR, RT-PCR and restriction digests.

Available in a pack of 1 with the mass of 50mg/500ml.

CAT No.: SMGBM-1001RL

Mass: 50mg/500ml

Application:

- Determine DNA fragments generated by PCR, RT-PCR and restriction digests

Features:

- Reduces set-up time
- Ready to use even if it is still frozen

Storage: Store at 2 to 8°C

Precaution:

- Once its thawed, store it at 2 to 8°C
- Avoid repeated freeze/thaw cycle